

# Targeted DNA/RNA Capture Magnetic Particles-SA PRODUCT DATA SHEET

# **Targeted DNA/RNA Capture Magnetic Particles-SA**

# Description

The streptavidin-biotin system has a very high binding affinity (Kd=10<sup>-15</sup>) and has been widely used in the biological field. This series of magnetic microspheres is a monodisperse microsphere product, which has the advantages of fast magnetic response, good suspension, large specific surface area, good stability, and small difference between batches. Surface modified streptavidin can bind biotinylated antibodies, nucleic acids, proteins and other ligand molecules quickly and efficiently, which is an ideal choice for capturing target molecules.

For custom sizes, formulations or bulk quantities please contact our customer service department. Website: www.abvigen.com Phone: +1 929-202-3014 Email: info@abvigenus.com

Product number	NMP	NMP	NMP	NMP	NMP	NMP
	0300SB	0600SB	1001SB	1002SB	1003SB	1005SB
Mean particle size	300	600 nm	1 μm	2 µm	3 µm	5 μm
	nm					
Free biological	/	> 800	> 700	1000	> 500	800
binding capacity						
(pmol/mg						
magnetic						
microspheres)						
Biotinized single	450	> 500	> 350	400	> 200	400
chain						
oligonucleotie						
binding amount						
(pmol/mg						
magnetic						
microspheres)						
Biotinized ds-DNA		≈ 20	10~20		5~10	
binding amount						

# Product information

1378 US-206 Ste 6-126, Skillman, NJ USA info@abvigenus.com Tel: 1-816-388- 0112 Fax: 1- 888-616-0161 Reserved Email:

© Abvigen Inc All Rights



(μg/mg magnetic						
microspheres)						
Biotinized IgG	15	> 20	≈ 15	20	≈ 10	15
binding amount						
(μg/mg magnetic						
microspheres)						
Magnetic	10 mg/mL					
microsphere						
concentration						
Preservation	10 mM Tris (pH	7.4), 0.1% BSA, 0	.05% Tween-20, 0	.1% ProClin 300		
solution						

# Characteristics

Concentration: 10 mg/ml

Particle Size: 0.3~5 um

Surface: Streptavidin

Preservative Solution: 10 mM Tris (pH 7.4), 0.1% BSA, 0.05% Tween-20, 0.1% ProClin 300

Storage Condition: Store at 2 - 8°C, do not freeze

Quality guarantee period: 24 months

# Scope of application

Biotin-modified nucleic acid probes are used for the capture and separation of nucleic acids. Binding biotin-modified antibodies or antigens for immunoassay, protein separation, etc.

#### Self-supplied material

#### **Recommended reagent**

DNA binding reagent	D-Binding Buffer: 10 mM Tris-HCl (pH 7.5), 1 mM EDTA, 2 M NaCl
	D-Washing Buffer: 5 mM Tris-HCl (pH 7.5), 0.5 mM EDTA, 1 M NaCl
RNA binding reagent	The following reagents should be prepared with Nuclease-free ddH <sub>2</sub> O:
	R-Binding Buffer: 0.1 M NaCl
	R-Washing Buffer: 0.1 M NaOH, 0.05 M NaCl
Protein binding reagent	PBST Buffer: PBS (pH 7.4), 0.05% Tween 20

1378 US-206 Ste 6-126, Skillman, NJ USA info@abvigenus.com Tel: 1-816-388- 0112 Fax: 1- 888-616-0161 Reserved Email:

© Abvigen Inc All Rights



	PBS/BSA Buffer: PBS (pH 7.4), 0.01% - 0.1% BSA
Nucleic acid release reagent	95% formamide solution (95% formamide + 5% 10 mM EDTA (pH 8.2))
Protein releasing reagent	0.1% SDS

#### Instruments and consumables

Vortex oscillators, vertical rotation mixers, magnetic racks, pipettes and matching heads, Nucleasefree eight row or suitable EP tubes.

#### **Experimental process**

The following experimental procedures are general procedures for streptavidin magnetic microspheres. Reagents and experimental procedures can be optimized according to requirements. According to the content of biotinylated molecules in the sample, the appropriate amount of magnetic microspheres is calculated by referring to the binding amount in the product information table.

#### Magnetic microsphere pretreatment

1. Balance the magnetic microspheres to room temperature, place them on the vortex oscillator for 30 s and mix them thoroughly, use a pipette to absorb the magnetic microspheres into the EP tube, place them on the magnetic rack, and discard the supernatant after the solution is clarified (about 1 min).

2. The EP tube was removed from the magnetic rack, and D-Washing Buffer (DNA capture) or R-Washing Buffer (RNA capture) or PBST Buffer (protein capture), which was twice the volume of the initial magnetic microsphere, was placed on the vortex oscillator and suspended for 10 s.

3. After a short centrifugation, place the EP tube on the magnetic rack and discard the supernatant after the solution is clarified.

4. DNA capture and protein capture are repeated (step 2,3) twice, and RNA capture is repeated once.

#### **DNA capture**

1. Add D-Binding Buffer (twice the volume of the original magnetic microspheres) to the pre-treated magnetic microspheres (discarded supernatant), and then add biotinized DNA samples equal to the volume of D-Binding Buffer. After fully mixing, Incubate at room temperature in a vertical rotating mixing machine for 15-30 min (mixing speed should not be too high, 20-30 rpm is recommended), after a short centrifugation, place the EP tube on a magnetic rack, and absorb and discard the supernatant after the solution is clarified.

1378 US-206 Ste 6-126, Skillman, NJ USA info@abvigenus.com Tel: 1-816-388- 0112 Fax: 1- 888-616-0161 Reserved



2. Add D-Washing Buffer twice the sample volume into EP tube, and blow 10 times with pipette to mix thoroughly. After a short centrifugation, the EP tube was placed on the magnetic rack, and the supernatant was sucked up after the solution was clarified.

3. Repeat Step 2.

4. According to the requirements of the follow-up experiment, add the appropriate buffer to resuspend the magnetic beads, and the magnetic beads can be used for subsequent operations.

#### **RNA** capture

1. An R-Binding Buffer of twice the original magnetic microsphere volume was added to the pretreated magnetic microsphere (discarded supernatant) and suspended on the vortex oscillator for about 10 s. After a short centrifugation, the EP tube was placed on the magnetic rack, and the supernatant was sucked up after the solution was clarified.

2. Add R-Binding Buffer twice the volume of the original magnetic microsphere, then add biotinized RNA sample equal to the volume of R-Binding Buffer, mix thoroughly, place on the magnetic rack at room temperature, and absorb the supernatant after the solution is clarified.

3. Add R-Binding Buffer solution, twice the volume of sample solution, into the EP tube, blow and mix with pipette 10 times, centrifuge briefly and place on the magnetic rack. After the solution is clarified, the supernatant is absorbed and discarded.

4. Repeat Step 3.

5. According to the requirements of the follow-up experiment, add the appropriate buffer to resuspend the magnetic beads, and the magnetic beads can be used for subsequent operations.

#### Protein capture

1. PBST Buffer of 2 times the volume of the original magnetic microspheres was added to the pretreated magnetic microspheres (discarded supernatant), and then the samples containing biotinated proteins to be captured were added. After being fully mixed, the samples were incubated in a vertical rotating mixing machine at room temperature for 15-30 min. After the solution is clarified, the supernatant is absorbed and discarded.

2. Add PBST Buffer or PBS/BSA Buffer, twice the original volume of magnetic microspheres, into the EP tube, blow 10 times with a pipette to mix thoroughly, after a short centrifuge, place the EP tube on the magnetic rack, and discard the supernatant after the solution is clarified.

3. Repeat Step 2.

1378 US-206 Ste 6-126, Skillman, NJ USA info@abvigenus.com Tel: 1-816-388- 0112 Fax: 1-888-616-0161 Reserved



4. Add PBST Buffer or PBS/BSA Buffer or other suitable buffer to re-suspension magnetic beads according to the requirements of subsequent experiments. Magnetic beads can be used for subsequent operations.

#### Biotin-modified nucleic acid/protein release

1. Separation of nucleic acid from streptavidin magnetic microspheres: After adding 95% formamide solution, heating at 65°C for 5 min or 90°C for 2 min, the nucleic acid is released.

2. Separation of protein from streptavidin magnetic microspheres: Add 0.1% SDS solution and boil for 5 min to separate protein.

#### Features

Superparamagnetic: Excellent resuspension
Hydrophilic surface: Low non-specific binding
Uniform diameter, CV < 5%: High reproducibility</li>
Coated with streptavidin (SA): Effectively binds biotinylated derivatives
Large scale production: Batch-to-batch consistency, superior quality with consistent test results

# Electron microscope image of Targeted DNA/RNA Capture Magnetic Particles-SA



#### Storage

#### This product should be stored at 2 - 8°C. DO NOT FREEZE.

1378 US-206 Ste 6-126, Skillman, NJ USA info@abvigenus.com Tel: 1-816-388- 0112 Fax: 1-888-616-0161 Reserved



#### Notes

1. Do not store the magnetic microsphere below 0°C, and balance the magnetic microsphere to room temperature before use (the volume is less than 10 mL for 30 min balance, greater than 10 mL, extend the balance time as appropriate).

2. Before absorbing magnetic microspheres, ensure that they are fully mixed to avoid bubbles caused by violent oscillation.

3. It is recommended to use a good quality pipette suction head and reaction tube to avoid losses caused by adhesion to magnetic microspheres and solution.

4. Avoid magnetic microspheres being centrifuged at high speed or placed on the magnetic rack for a long time.

5. Avoid magnetic microspheres in a liquid free state for a long time. If there are many samples, it is recommended to discard the supernatant in batches and add the re-suspension reagent.

6. If there is a phenomenon of magnetic microspheres sticking to the wall in a high-salt solution,0.01%-0.1% Tween 20 can be added to the solution in advance.

7. This product is for scientific research only.

# **Ordering Information**

Website: <u>www.abvigen.com</u> Phone: +1 929-202-3014 Email: <u>info@abvigenus.com</u>

1378 US-206 Ste 6-126, Skillman, NJ USA info@abvigenus.com Tel: 1-816-388- 0112 Fax: 1- 888-616-0161 Reserved