



Gold Nanoparticles Conjugation Kit-Oligo

PRODUCT DATA SHEET

Gold Nanoparticles Conjugation Kit-Oligo

Description

Oligonucleotides are nucleic acid polymers with the potential to treat or manage a wide range of diseases. Specifically, Oligonucleotides can be used to modulate gene expression via a range of processes including RNAi, target degradation by RNase H-mediated cleavage, splicing modulation, non-coding RNA inhibition, gene activation and programmed gene editing. As such, these molecules have potential therapeutic applications for myriad indications. In addition, it is worth noting that the highly specific molecular recognition properties of oligonucleotides enable them to combine with the unique optical properties of gold nanoparticles, thereby developing biosensors for application in the field of biotechnology.

Abvigen Gold Nanoparticles Conjugation Kit-Oligo have been optimized for high efficiency one-step conjugation of thiolated oligonucleotides directly to the gold surface of particles with diameters up to 100nm. The kit contains ready-to-use pre-made mixtures. No activation, manipulation, or time consuming "salt-aging" steps are required for conjugation. Simply mix your reduced thiol-modified oligonucleotide with the pre-activated gold nanoparticles supplied with the kit. Conjugation of the oligonucleotide is achieved by the formation of a strong and stable gold-thiol bond without any additional linkers.

Kits are available in convenient 3 or 10 small-scale reaction formats allowing multiple to be conjugated simultaneously. The precisely engineered gold surface on our Oligo gold nanoparticles results in high conjugation efficiency and stable conjugates while minimizing non-specific binding in your assay.

For custom sizes, formulations or bulk quantities please contact our customer service department.

Website: www.abvigen.com **Phone:** +1 929-202-3014 **Email:** info@abvigenus.com

Characteristics

Gold surface: Oligo

Core diameter: Available with diameters from 5 nm ~ 100 nm

Optical density (OD): OD=2 when the contents of each vial is dissolved to a final volume of 1 ml.



Features & Benefits

Allows conjugation of oligonucleotides to gold nanoparticles with sizes between 5 nm ~ 100 nm.

Fast and convenient one-step conjugation reaction with no pre-activation requirements or manipulation of the gold nanoparticles.

No time-consuming "salt-aging" procedures.

Results in a thiol-oligo conjugated directly to the gold surface without any linkers.

Optimized for use in crosslinking based lateral flow applications.

Procedure

Reduction of thiol-modified oligonucleotides (e.g. trityl-S-S-Oligo)

1. Prepare a 0.15 M sodium phosphate buffer, pH 8.5 supplemented with 0.1 M DTT.

Note: pH is important for proper reduction of oligonucleotide.

2. Dissolve lyophilized oligonucleotide to a final concentration of 500 μ M in H₂O.

3. Mix 50 μ l of dissolved oligonucleotide with 450 μ l sodium phosphate buffer.

4. Incubate 1-2 h at room temperature to reduce oligonucleotide.

5. Separate reduced oligonucleotide from trityl-SH and DTT using a NAP 5 column operated in H₂O, GE Healthcare.

6. Final eluate from NAP 5 column will be 1 ml in H₂O with an approximate concentration of 25 μ M.

Note: Exact concentration of final eluate can be measured with UV-VIS spectroscopy by measuring the absorbance at 260 nm.

Conjugation of thiolated oligonucleotide to Oligo gold nanoparticles

1. Resuspend one vial of lyophilized Oligo gold nanoparticle with 740 μ l of H₂O.

2. Transfer into a 1.5 ml microcentrifuge tube.

3. Add 160 μ l of reduced thiolated oligonucleotide at 7.5 μ M (0.0075 nmol/ μ l)* in H₂O as prepared above and incubate for at least 1 h at room temperature.

*Note: 7.5 μ M oligonucleotide is a good starting concentration, but if aggregation or poor sensitivity is observed, the following oligonucleotide concentrations can be attempted for a given particle size range (based on a 30nt oligonucleotide):

Particle size (nm)	5	10-20	30-100
[oligonucleotide] (μ M)	5-50	5-25	5-15

4. Add 100 μ l of 1 M NaCl.



5. Incubate for at least 1 h at room temperature to allow binding of the oligonucleotide to the gold surface.

Note: Longer incubation times may improve surface coverage.

6. Centrifuge at the appropriate speed for your particular gold nanoparticle size (see table I) for 30 min to pellet your oligonucleotide gold conjugate.

7. Remove supernatant.

8. Resuspend conjugate in 200 µl of storage buffer.

The optical density of the particles should be 10 if a 100% recovery has been achieved.

Common storage buffer: 10 mM sodium phosphate buffer, pH 7.0, 100 mM NaCl and 0.01% (w/v) NaN₃.

9. Measure optical density with a spectrophotometer and adjust concentration as desired.

10. Store conjugate at +4°C.

Table I. Appropriate G forces for centrifugation of gold nanoparticles. Note that recommended conditions are for a volume of 1ml and centrifugation using a microcentrifuge, except for 5 nm gold nanoparticles that require an ultracentrifuge.

Size (nm)	Speed (g)	Time (min)
5	100,000	30
10	17,000	60 (~50% recovery)
15	17,000	30
20	6,500	30
30	4,500	30
40	2,500	30
50	2,000	30
60	1,125	30
80	600	30
100	400	30

Storage

Store at -20°C. Stable for at least 3 months if stored as specified.

**Note**

This product is for R&D use only, not for drug, household, or other uses.

Ordering Information

Website: www.abvigen.com

Phone: +1 929-202-3014

Email: info@abvigenus.com