

CD8 Magnetic Nanoparticles

Description

CD8 magnetic nanoparticles are homogeneous and biocompatible superparamagnetic nanoparticles coated with highly specific monoclonal antibodies, which can be used for the isolation of human CD8⁺T cells. CD8 magnetic nanoparticles are small in size, good in suspension, and specifically target CD8 positive cells. Cells bound with magnetic beads can be enriched and purified by magnetic separation, and then cells of high purity can be obtained by positive sorting or removal of CD8 cells. The separation can be completed by manual or automatic sorting equipment. Obtain high purity, high activity and good functioning human CD8⁺T cells.

Abvigen offers high quality CD8 magnetic nanoparticles. The product has high repeatability between batches, which can meet the needs of various customers for personalized materials such as research and development, testing and production.

For custom sizes, formulations or bulk quantities please contact our customer service department. Website: <u>www.abvigen.com</u> Phone: +1 929-202-3014 Email: <u>info@abvigenus.com</u>

Characteristics

Type: CD8 Magnetic Nanoparticles Surface group: Anti-CD8 Dispersing solvent: Preservative solution Size: 2 mL Storage condition: Store at 2-8°C away from light to avoid freezing.

Operation Case

Material preparation

- 1. Magnetic sorting column and magnetic field (magnetic sorting rack);
- 2. Separation buffer: 0.01 M PBS+0.5% BSA+ 2 mM EDTA, pH=7.2, refrigerated at 2-8°C for use;
- 3. Gently mix CD8 nanobeads before use.

Cell preparation

1378 US-206 Ste 6-126, Skillman, NJ USA Tel: 1-816-388- 0112 Fax: 1- 888-616-0161



1. When human peripheral blood is used as the isolation sample, peripheral blood mononuclear cells (PBMC) should be isolated using lymphocyte isolation solution;

2. Resuspension PBMC with separation buffer to 1×10⁷ PCS /mL, refrigerated at 2-8°C for use.

Magnetic separation

Take a single isolation of 1×10^7 cells as an example:

1. Centrifuge the prepared cell suspension at 300 g for 10 min and discard the supernatant;

2. The cells were precipitated and suspended in 80 μ L separation buffer, and 20 μ L CD8 nanobeads were added, thoroughly mixed, and incubated at 2-8°C for 15 min.

3. After incubation, add 1 mL separation buffer, centrifuge 300 g for 10 min, discard the supernatant, and re-suspend the cells with 1 mL separation buffer;

4. Place the magnetic sorting column in the magnetic field, add 1 mL separation buffer to wash, and add the cell suspension incubated by CD8 nanoparticle into the magnetic sorting column;

5. When the liquid in the magnetic sorting column flows out, the separation buffer is added twice, adding 1 mL each time, and all the effluent is collected. The effluent is the PBMC that removes CD8⁺T cells;

6. Remove the magnetic sorting column from the magnetic field, place it on a suitable collection tube, add 1 mL separation buffer to the magnetic sorting column, and use a push rod to flush out the cells trapped in the sorting column, namely CD8⁺T cells.

Note:

(1) The number of cells isolated in a single time is recommended to be no more than 1×10^7 , and when a higher number of cells is used, the amount of reagents and the total volume of the reaction should be increased accordingly;

(2) Keeping the cells operating at 2-8°C can effectively reduce non-specific binding;

(3) The purity of target cells can be improved by using the new magnetic sorting column.

Advantages

Small size Good suspension

Applications

The isolation of human CD8⁺T cells



Storage

Store at 2-8°C away from light to avoid freezing.

Ordering Information

Website: <u>www.abvigen.com</u> Phone: +1 929-202-3014 Email: <u>info@abvigenus.com</u>