



PRODUCT DATA SHEET

Blue Fluorescent SiO₂ Magnetic Particles-Biotin

Description

Fluorescent Particles are prepared by either incorporation of selected fluorophores into monodisperse polystyrene particles by means of swelling processes or by copolymerization of styrene with various organic fluorescent dyes, these processes generate fluorophores labeled polystyrene particles with favorable properties. As a result, a wide variety of fluorescent particles can be prepared ranging in size, type of fluorophore, fluorescence intensity and surface functional groups. Most of the fluorophores chosen for use in the preparation of Fluorescent Particles are water insoluble and therefore are very stable. These fluorophores, once incorporated into the particles, do not leach and their color and fluorescence remains stable for long periods of time under proper storage conditions. Fluorescent particles with single or multiple fluorophores are available in various sizes, emission spectra and combinations. Many are suitable for uses in flow cytometry, fluorescence microscopy, phagocytosis studies and cell labeling.

For custom sizes, formulations or bulk quantities please contact our customer service department.

website: www.abvigen.com Phone: +1 929-202-3014 Email: info@abvigenus.com

Characteristics

Concentration: 10 mg/ml

Size: 10 ml; 20 ml

Surface: Biotin

Shape: Spherical

Composition: Magnetic SiO₂ Particles

Standard deviation: CV<5%

Excitation: 360-435 nm

Emission: 450 nm

Buffer: PBS

Form: Suspension

Store: Storage at 2 - 8 °C



Storage

This product should be stored at 4°C. **DO NOT FREEZE.**

Ordering Information

website: www.abvigen.com

Phone: +1 929-202-3014

Email: info@abvigenus.com