



Time Resolved Fluorescent Particles-Biotin

PRODUCT DATA SHEET

Time Resolved Fluorescent Particles-Biotin

For custom sizes, formulations or bulk quantities please contact our customer service department.

website: www.abvigen.com Phone: +1 929-202-3014 Email: info@abvigenus.com

Characteristics

Core diameter: 200 nm; 300 nm; 400 nm

Size: 10 ml; 20 ml; 50 ml; 100 ml

Concentration: 10 mg/ml

Surface: Biotin (150-200 $\mu\text{mol/g}$)

PDI: ≤ 1.05

Shape: Spherical

Density: 1.03 g/ccm

Composition: Europium (Eu) Time Resolved Fluorescent Particles

Excitation: 360 nm

Emission: 615 nm

Buffer: DI Water

Form: Suspension

Storage

This product should be stored at 4°C. **DO NOT FREEZE.**

Ordering Information

website: www.abvigen.com

Phone: +1 929-202-3014

Email: info@abvigenus.com

1378 US-206 Ste 6-126, Skillman, NJ USA

info@abvigenus.com

Tel: 1-816-388- 0112 Fax: 1- 888-616-0161

Reserved

Email:

© Abvigen Inc All Rights