



## Basic mRNA Purification Kit PRODUCT DATA SHEET

### Basic mRNA Purification Kit

**Cat No: AKIT-Basics**

#### Description

The base mRNA purity test kit is a prefabricated buffer specially designed and developed based on Oligo-dT coated magnetic beads to provide the best stable conditions for mRNA purification. It is used for rapid purification of high purity and complete mRNA from total RNA extracted from animal tissues, plants or cells, which can be used for various downstream experiments in molecular biology.

For custom sizes, formulations or bulk quantities please contact our customer service department.

**Website:** [www.abvigen.com](http://www.abvigen.com) **Phone:** +1 929-202-3014 **Email:** [info@abvigenus.com](mailto:info@abvigenus.com)

#### Kit Composition

Binding buffer	5 mL
Washing solution	5 mL
RNase-Free H <sub>2</sub> O	5 mL

#### Operation Process

This kit is only used with Oligo-dT coated magnetic beads. For details, see our Oligo-dT coated magnetic beads or mRNA purification kit instructions.

#### Ordering Information

**Website:** [www.abvigen.com](http://www.abvigen.com)

**Phone:** +1 929-202-3014

**Email:** [info@abvigenus.com](mailto:info@abvigenus.com)